


MC3P 612-Y



on a scale of 5:1

Technical Information

Type	Posterior Chamber IOL with Blue Light Protection (Three-Piece, Foldable, 360° Sharp Posterior Optic Edge)
Optic Material	Hydrophilic Acrylic with UV-Absorber and Yellow Colouring Agent
Water Content	26% at 35°C
Haptic Material	Highly Flexible Colourized PES
Haptic Shape	C-Loop
Optic Shape	Aspheric Posterior Surface (Aberration Free)*
	Biconvex
Dioptr Range	10.0 to 30.0 D
Dioptr Steps	0.5 D
A-Constant (est.)**	118.6

Dr. Schmidt Intraocularlinsen GmbH, Germany

Westerwaldstr. 11-13
53757 Sankt Augustin

Phone: +49 (0) 2241 / 25 787-0
Fax: +49 (0) 2241 / 25 787-88

mail@dr-schmidt-iol.de
www.dr-schmidt-iol.de

* The word „aberration“ as used in this document is interpreted as „spherical aberration“

** For current status and more details please refer to www.dr-schmidt-iol.de/a-constants.pdf