


MC 5812 AS



on a scale of 5:1

Technical Information

<b>Type</b>	Posterior Chamber IOL (Single-Piece, Foldable)		
<b>Material</b>	Hydrophilic Acrylic with UV-Absorber		
<b>Water Content</b>	26% at 35°C		
<b>Haptic Shape</b>	C-Loop		
<b>Optic Shape</b>	Concave-convex	Biconvex	
<b>Diopter Range</b>	0.0 to 9.0 D	10.0 to 30.0 D	31.0 to 35.0 D
<b>Diopter Steps</b>	1.0 D	0.5 D	1.0 D
<b>A-Constant (est.)*</b>	118.4		

Dr. Schmidt Intraocularlinsen GmbH, Germany

Westerwaldstr. 11-13  
53757 Sankt Augustin

Phone: +49 (0) 2241 / 25 787-0  
Fax: +49 (0) 2241 / 25 787-88

mail@dr-schmidt-iol.de  
www.dr-schmidt-iol.de

\* For current status and more details please refer to [www.dr-schmidt-iol.de/a-constants.pdf](http://www.dr-schmidt-iol.de/a-constants.pdf)