


MC 6105



on a scale of 5:1

Technical Information

Type	Posterior Chamber IOL (Single-Piece, Foldable)	
Material	Hydrophilic Acrylic with UV-Absorber	
Water Content	26% at 35°C	
Optic	Aspheric Anterior Surface (Aberration Free)* Posterior Surface with 360° Lens Epithelial Cell Barrier	
Haptic Shape	4-Loop	
Optic Shape	Concave-convex	Biconvex
Diopter Range	0.0 to 9.0 D	10.0 to 30.0 D
Diopter Steps	1.0 D	0.5 D
A-Constant (est.)**	118.7	

Dr. Schmidt Intraocularlinsen GmbH, Germany

Westerwaldstr. 11-13
53757 Sankt Augustin

Phone: +49 (0) 2241 / 25 787-0
Fax: +49 (0) 2241 / 25 787-88

mail@dr-schmidt-iol.de
www.dr-schmidt-iol.de

* The word „aberration“ as used in this document is interpreted as „spherical aberration“

** For current status and more details please refer to www.dr-schmidt-iol.de/a-constants.pdf