


MC 6125 AS



on a scale of 5:1

Technical Information

<b>Type</b>	Posterior Chamber IOL (Single-Piece, Foldable)		
<b>Material</b>	Hydrophilic Acrylic with UV-Absorber		
<b>Water Content</b>	26% at 35°C		
<b>Optic</b>	Aspheric Anterior Surface (Aberration Free)* Posterior Surface with 360° Lens Epithelial Cell Barrier		
<b>Haptic Shape</b>	C-Loop		
<b>Optic Shape</b>	Concave-convex	Biconvex	
<b>Diopter Range</b>	-10.0 to 9.0 D	10.0 to 30.0 D	31.0 to 50.0 D
<b>Diopter Steps</b>	1.0 D	0.5 D	1.0 D
	<b>More details on diopter categories see page 2</b>		
<b>A-Constant (est.)**</b>	118.4		

Dr. Schmidt Intraocularlinsen GmbH, Germany

Westerwaldstr. 11-13  
53757 Sankt Augustin

Phone: +49 (0) 2241 / 25 787-0  
Fax: +49 (0) 2241 / 25 787-88

mail@dr-schmidt-iol.de  
www.dr-schmidt-iol.de

03/2011

## Dioptr Category – MC 6125 AS

	Standard	0.0 to 9.0 D 10.0 to 30.0 D	in 1.0 D steps in 0.5 D steps
	Special	-6.0 to -1.0 D 31.0 to 40.0 D	in 1.0 D steps in 1.0 D steps
	Customized	-10.0 to -7.00 D 41.0 to 50.0 D	in 1.0 D steps in 1.0 D steps
	Customized plus	> 50.0 D	on request

The dioptr category does not indicate any conclusion related to availability. The delivery time of single powers might vary and is subject to change.