


MS 614 ASP-Y



Technical Information

Type	Posterior Chamber IOL with Blue Light Protection (Three-Piece, Foldable) for Sulcus Fixation
Optic Material	Silicone Elastomer with UV-Absorber and Yellow Colouring Agent
Haptic Material	High Molecular PMMA
Haptic Shape	Modified C-Loop
Optic Shape	Spheric Anterior Surface, Aspheric Posterior Surface (Aberration Correcting)*
	Biconvex
Diopter Range	7.0 to 31.0 D
Diopter Steps	0.5 D
A-Constant (est.)**	118.0

Dr. Schmidt Intraocularlinsen GmbH, Germany

Westerwaldstr. 11-13
53757 Sankt Augustin

Phone: +49 (0) 2241 / 25 787-0
Fax: +49 (0) 2241 / 25 787-88

mail@dr-schmidt-iol.de
www.dr-schmidt-iol.de

* The word „aberration“ as used in this document is interpreted as „spherical aberration“

** For current status and more details please refer to www.dr-schmidt-iol.de/a-constants.pdf