


MS 614 T



on a scale of 5:1

Technical Information

Type	Toric Posterior Chamber IOL (Three-Piece, Foldable) for Sulcus Fixation
Optic Material	Silicone Elastomer with UV-Absorber
Haptic Material	High Molecular PMMA
Haptic Shape	Modified C-Loop
Optic Shape	Spheric Anterior Surface, Toric Posterior Surface
Sphere Range	in 0.5 D steps: 15.0 to 25.0 D in 1.0 D steps: -30.0 to 14.0 D / 26.0 to 34.0 D
Cylinder Range	in 0.5 D steps: 1.0 to 3.0 D in 1.0 D steps: 4.0 to max. 30.0 D*
	Availability depending on combination of sphere and cylinder. More details on dioptric categories see page 2
A-Constant (est.)*	118.0
Cylinder Axis	The Two Thin Diametrically Opposed Lines in the Periphery of the Optic Indicate the Axis of the Plus-Cylinder


Dr. Schmidt Intraocularlinsen GmbH, Germany

Westerwaldstr. 11-13
53757 Sankt Augustin

Phone: +49 (0) 2241 / 25 787-0
Fax: +49 (0) 2241 / 25 787-88

mail@dr-schmidt-iol.de
www.dr-schmidt-iol.de

Dioptr Category – MS 614 T



The dioptr category does not indicate any conclusion related to availability. The delivery time of single powers might vary and is subject to change.