


MS 714 PB Diff-Y



on a scale of 5:1

Technical Information

Type	Multifocal Posterior Chamber Add-On IOL with Blue Light Protection (Three Piece, Foldable) for Sulcus Fixation in Pseudophakic Eyes
Optic Material	Silicone Elastomer with UV-Absorber and Yellow Colouring Agent
Haptic Material	High Molecular PMMA
Haptic Shape	Modified C-Loop
Optic Shape	Diffractive Anterior Surface, Spheric Posterior Surface (≈ Aberration Free)*
	Convex-concave with Round Anterior Optic Edge to Prevent Iris Irritation
Dioptr Range	-6.0 to 6.0 D
Dioptr Steps	0.5 D
	More details on dioptr categories see page 2
Near Addition	+ 3.5 D (on IOL Plane)
A-Constant (est.)	Not Applicable (by Application Principle)

Dr. Schmidt Intraocularlinsen GmbH, Germany

Westerwaldstr. 11-13
53757 Sankt Augustin

Phone: +49 (0) 2241 / 25 787-0
Fax: +49 (0) 2241 / 25 787-88

mail@dr-schmidt-iol.de
www.dr-schmidt-iol.de

Dioptr Category – MS 714 PB Diff-Y

	Standard	0.0 D	
	Special	-6.0 to -0.5 D 0.5 to 6.0 D	in 0.5 D steps in 0.5 D steps

The dioptr category does not indicate any conclusion related to availability. The delivery time of single powers might vary and is subject to change.