


MS 714 PB Diff



Technical Information

<b>Type</b>	Multifocal Posterior Chamber Add-On IOL (Three Piece, Foldable) for Sulcus Fixation in Pseudophakic Eyes
<b>Optic Material</b>	Silicone Elastomer with UV-Absorber
<b>Haptic Material</b>	High Molecular PMMA
<b>Haptic Shape</b>	Modified C-Loop
<b>Optic Shape</b>	Diffractive Anterior Surface, Spheric Posterior Surface ( $\approx$ Aberration Free)*
	Convex-concave with Round Anterior Optic Edge to Prevent Iris Irritation
<b>Dioptr Range</b>	-6.0 to 6.0 D
<b>Dioptr Steps</b>	0.5 D
	<b>More details on dioptr categories see page 2</b>
<b>Near Addition</b>	+ 3.5 D (on IOL Plane)
<b>A-Constant (est.)</b>	Not Applicable (by Application Principle)

Dr. Schmidt Intraocularlinsen GmbH, Germany

Westerwaldstr. 11-13  
53757 Sankt Augustin

Phone: +49 (0) 2241 / 25 787-0  
Fax: +49 (0) 2241 / 25 787-88

mail@dr-schmidt-iol.de  
www.dr-schmidt-iol.de

## Dioptr Category – MS 714 PB Diff

	Standard	0.0 D	
	Special	-6.0 to -0.5 D 0.5 to 6.0 D	in 0.5 D steps in 0.5 D steps

The dioptr category does not indicate any conclusion related to availability. The delivery time of single powers might vary and is subject to change.