


MS 714 PB-Y



Technical Information

Type	Posterior Chamber Add-On IOL with Blue Light Protection (Three-Piece, Foldable) for Sulcus Fixation in Pseudophakic Eyes
Optic Material	Silicone Elastomer with UV-Absorber and Yellow Colouring Agent
Haptic Material	High Molecular PMMA
Haptic Shape	Modified C-Loop
Optic Shape	Convex-concave with Round Anterior Optic Edge to Prevent Iris Irritation (≈ Aberration Free)*
Dioptr Range	-6.0 to 6.0 D
Dioptr Steps	0.5 D
A-Constant (est.)	Not Applicable (by Application Principle)

Dr. Schmidt Intraocularlinsen GmbH, Germany

Westerwaldstr. 11-13
53757 Sankt Augustin

Phone: +49 (0) 2241 / 25 787-0
Fax: +49 (0) 2241 / 25 787-88

mail@dr-schmidt-iol.de
www.dr-schmidt-iol.de

* The word „aberration“ as used in this document is interpreted as „spherical aberration“